

## **Product Information**

# Cellvento<sup>™</sup> Feed-200 Chemically defined cell culture feed

## **Product description**

Cellvento<sup>™</sup> Feed-200 is a chemically defined feed formulation containing only components of non-animal origin. The product is intended for use in the development and manufacturing of biotherapeutics in Chinese Hamster Ovary (CHO) cell-based expression systems.

Cellvento™ Feed-200 contains concentrated amino acids. vitamins, salts, and trace elements and is supplied in dry powder format.

#### **Application**

Cellvento<sup>™</sup> Feed-200 has been designed to support optimal cell growth and performance of CHO-S cells, but may be suitable for use with other CHO cell lines. It is intended for use as a companion feed supplement with the production medium Cellvento™ CHO-200 in fed-batch manufacturing processes.

This product is intended for research or further manufacturing but not for human or therapeutic use.

## Reconstitution method to generate 3 L Cellvento<sup>™</sup> Feed-200

- 1. Slowly add 347 grams of powder to 2.5 L of Milli-Q® or similar cell culture grade water in an appropriately sized container. Rinse feed container as necessary to remove remaining powder.
- 2. Vigorously mix for 45-60 minutes until fully dissolved.
- 3. Add cell culture grade water to reach a final volume of 3 L. Confirm final pH of 5.3-5.9.
- 4. Measure the osmolality of the solution. Final osmolality should be 1150 +/-50 mOsmol/kg. Sterilize by membrane filtration using a 0.22 μm Millipore Express® PLUS or Durapore® membrane filter (bottle cap or capsule filter).
- 5. Store at 2-8 °C protected from light. Reconstituted liquid Cellvento™ Feed-200 is stable for 30 days. When a bottle is opened, liquid feed is stable for max. 2 weeks.

#### **Storage**

Dry powder should be stored at 2-8 °C protected from light. Do not use after expiration date.

#### Shelf life

12 months



#### Ordering information for companion Cellvento™ Feed-200

| Catalog number | Product name   | Pkg. size | Equivalent |
|----------------|--|-----------|------------|
| 1.01883.0003   | Cellvento™ Feed-200<br>Chemically defined cell<br>culture feed | 0.347 kg  | 3 liters   |
| 1.01883.0010   | Cellvento™ Feed-200<br>Chemically defined cell<br>culture feed | 1.158 kg  | 10 liters  |
| 1.01883.0050   | Cellvento™ Feed-200<br>Chemically defined cell<br>culture feed | 5.791 kg  | 50 liters  |

### Ordering information for Cellvento™ CHO-200 medium

| Catalog number | Product name  | Pkg. size | Equivalent |
|----------------|---|-----------|------------|
| 1.01885.0010   | Cellvento™ CHO-200<br>Chemically defined cell<br>culture medium | 0.232 kg  | 10 liters  |
| 1.01885.0100   | Cellvento™ CHO-200<br>Chemically defined cell<br>culture medium | 2.325 kg  | 100 liters |

#### Ordering information for cell culture additives

| Catalog number       | Product name   | Pkg. size |
|----------------------|--|-----------|
| 1.00286.1000         | L-Glutamine suitable for use as excipient EMPROVE® exp DAB, USP                                      | 1 kg      |
| 1.37013.2500         | Sodium hydrogen carbonate suitable for biopharmaceutical production EMPROVE® bio Ph Eur, BP, USP, JP | 2.5 kg    |
| Available on request | HT Supplement (50 ×)   | 100 mL    |
| 1.02735.0100         | L-Cysteine hydrochloride monohydrate<br>suitable for use as excipient EMPROVE®<br>exp Ph Eur, USP    | 0.1 kg    |
| 1.02413.0100         | L-Tyrosine disodium salt dihydrate for cell culture media  | 0.1 kg    |
| 1.02415.0400         | D(+) Glucose anhydrous for cell culture media  | 0.4 kg    |

#### Ordering information for aseptic filters

| Catalog number | Product name  | Pkg. size |
|----------------|---|-----------|
| GPWP02500      | Millipore Express® PLUS Membrane, 0.22<br>μm, 25 mm | 100       |
| GVWP02500      | Durapore® Membrane, 0.22 μm, 25 mm                  | 100       |

To find out more about Cellvento™ CHO media platform products, visit www.merckmillipore.com/cellvento

## For more information and documentation please contact:

Phone: +49 6151 72-0

Email: pcs.salessupportEU@merckgroup.com



Merck Millipore Merck KGaA Frankfurter Str. 250 64293 Darmstadt, Germany

## www.merckmillipore.com

Merck Millipore, the M mark, Durapore, EMPROVE, Millipore Express and Milli-Q are registered trademarks of Merck KGaA, Darmstadt, Germany. Cellvento is a trademark of Merck KGaA, Darmstadt, Germany. © 2015 Merck KGaA, Darmstadt, Germany. All rights reserved.

The typical technical data above serve to generally characterize the cell culture media in industry-relevant expression systems. The product information is available separately from the website <a href="https://www.merckmillipore.com">www.merckmillipore.com</a>

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.